


STARSYSTEM

Adept Medical introduces the STARSystem. A combination of the STARBoard, STARSupport and STARTable to become a complete, clinically engineered solution for radial access procedures.

STARSystem provides unequalled patient comfort and ease of use for clinicians. Designed to work with any cardiac or radiology procedure table, the STARSystem is quick to install and adjust to variations in patient size. Set-up is achieved in seconds with no consumable components required.



reddot award 2015
honourable mention



KEY FEATURES

Clinician design

Scatter Shielding


Fits all tables

Quick set-up

Easy positioning

Carbon fibre

Ergonomic surface



STARBOARD

The key function of the STARBoard is its ability to present the patients wrist in the extended position whilst access is achieved, then simply return it to a more relaxed, medially-rotated position, with the arm at their side for the duration of the procedure. Adjustments can be made at any point during the procedure beneath the sterile drapes.

A simple support strap around the upper forearm eliminates the need for taping restless patients.

Crafted in carbon fibre for superior strength, radiolucency and durability, the STARBoard is extremely light weight and compact.

The unique, one-piece, "foldaway" design makes it easy for nursing staff to handle and takes up minimal storage space.

STARSUPPORT

STARSupport connects to the STARBoard after radial access is gained facilitating left arm procedures for superior patient comfort.





STARTABLE

The STARTable provides clinicians with an adjustable work surface and protection from scatter radiation. For a full report please visit www.adeptmedical.com/starsystem

With STARTable

Without STARTable



SPECIFICATIONS

STARSystem
Product code: AM0160C

STARBoard + STARSupport
Product code: AM0100D

STARTable
Product code: AM0140C

STARBoard weight: 0.8Kg

STARSupport weight: 0.2Kg

STARTable weight: 2.3Kg

EC REP



EC REPRESENTATIVE

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

For more information please visit
www.adeptmedical.com/starsystem