


*"The Thompson Hip System has an elegant design that optimizes the synergy between instrumentation and surgical technique allowing **unimpeded visualization** while performing muscle sparing (or minimally invasive) total hip arthroplasty."*

- Mark J. Powers, MD, FACS, FAAOS



Uncompromised Exposure™ for Hip with the Thompson Retractor

Table-Mounted Hip Retractors

Firm Hold with Fewer Hands Required

Thompson table-mounted hip retractors are ideal for orthopedic surgeons performing total hip replacement. These specialized systems utilize rigid handles or flexible tethers to attach retractor blades to a stable, table mounted frame which eliminates the need for extra staff to hold retractors. A sterile mount, low profile design, and streamlined set up saves time and space.

Thompson Retractor Exclusive Features: Thompson is dedicated to offering the most forward thinking, proven, and high quality features. These exclusive, often patented, features define Thompson and allow for Uncompromised Exposure.

“By improving the direct visualization of the native landmarks during Direct Anterior Approach THA, the Thompson Hip Retractor System provides the surgeon with improved confidence for consistent, accurate component placement.”

- Lee E. Rubin, MD



Anatomically Designed Frame

Low profile, independent frame arms contour to patient anatomy.



Compatibility

Thompson Hip Retractor systems are compatible with most hip retractor blades.



Sterile Mount Rail Clamp

Thompson's table mounted frame eliminates retractor shifting. Top control allows attachment to table, as well as adjustments, to be achieved within the sterile field.



Universal Blade Adapter Handles

Quick set up, with rigid attachment along flat blade surface.



Tethers

Fast, easy set up and flexible attachment.

Hip Retractors: Attachment Options

Thompson offers several retractor blade attachment options, from quick flexible tethers, to a variety of rigid handles.

Tethers

Flexible tethers attach quickly and securely while remaining flexible to allow slight movements when needed.

#42105 (5-pack)



Universal Blade Adapter Handle

Adapter Handles attach quickly along flat blade surface, providing a secure, rigid hold.

#42126BP



Orthopedic Clip-On Hook Handle

Fast, easy option for blade attachment and adjustments. Compatible with most hip blades containing holes.

#42126OHP



Orthopedic Handle and Knob

Firmly attach to most orthopedic blades containing holes, with the use of a universal attachment knob.

#42126OP (Handle)

#60025 (Knob)

Hip Retractors: Blades

Choose from Thompson's selection of retractor blades, or use your own retractors if preferred.



Curved Hohmann Wide 12"

#45307



Femoral Neck Elevator

#45300



Curved Hohmann Narrow 12"

(Color coded for easy identification)

#45306C



Hohmann 11"

#45305



Curved Dual Prong 12"

#45308



Dual Prong 11"

#45309



Curved Cobra 12 1/2"

(Color coded for easy identification)

#45302C



Cobra 11"

#45301

Thompson Retractor

Hip Retractor Add-Ons

Complete Your System

In addition to our Hip Retractor Frames, Thompson offers a variety of products to facilitate retraction during hip procedures.


Femur Elevator

Provides a fast connection and controlled, precise elevation using the micro-adjustable elevator.



S-Lock® Articulating Arm + Long Hibbs Blade

May be used for initial access and closing as desired.



Muscle Retractor Add-On Kit

Provides stable muscle retraction using Micro-Adjustable Handles and a variety of blade depths. Especially helpful for posterior approach.



Rail Extender

Allows for compatibility with Hana® Table and leg positioning systems.




Clip-On Cleats


Allow for attachment of tethers onto smooth frame arms where needed, with easy on and off "pinch" joint.

Hip Retractor Suggested Set Ups


Tethers Set Ups



Anterior Approach






Lateral Approach



Posterior Approach

Handles Set Ups



For a Free Trial Call Today*

1.800.227.7543

Rev I
053019
mg hipsg



Uncompromised Exposure™

VISIT US ONLINE: **THOMPSONSURGICAL.COM**

CONTACT

Address

Thompson Surgical Instruments, Inc.
10341 E Cherry Bend Rd.
Traverse City, MI 49684
USA

Phone + Fax

Toll Free 1-800-227-7543
Phone 231-922-0177
Fax 231-922-0174

Online

Email info@thompsonsurgical.com
Website www.thompsonsurgical.com

Social



Thompson
SURGICAL INSTRUMENTS

**Made In
USA**

* Free trial valid for U.S. customers only. Customers outside the U.S.A., please call +1-231-922-0177 for availability.

© 2019 Thompson Surgical Instruments, Inc. Traverse City, Michigan. Printed in U.S.A.

© S-Lock®, PLA®, and the "T Circle" logomark are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064. Other patents pending.

Note: As we continually strive to provide the best products possible, some of the images in this document may appear slightly different from the product received.